

Motion Custom

User Guide

SIEMENS

Content

Your Motion Custom hearing

instrument	4
Intended use	6
General notes on safety	7
Handling batteries	9
Low battery indicator	9
Inserting the battery	10
Inserting and removing the hearing instrument	11
Inserting	11
Removing	11
Turning your hearing instrument on and off	12
Via battery compartment	12
Via program button	12
Via volume control with on/off-function	12
Via remote control	12
Power-on delay	13
Changing the hearing program	14

Via program button	14
Via remote control	14
Changing the volume	15
Automatic	15
Via volume control	15
Via remote control	15
Functionalities	16
e2e wireless	16
Maintenance and care	19
Wax protection	20
Accessories	21
Tek wireless enhancement	21
Remote control	21
Troubleshooting	22
Technical information	23
Disposal information	25

Your Motion Custom hearing instrument




IT / ITE

(battery type: 13 or 312*)



CT / ITC / HS

(battery type: 312 or 10*)



- ① Battery compartment
- ② Microphone
- ③ Second Microphone
TwinMic™ (optional,
not for CIC)
- ④ Program button (optional)
- ⑤ Volume control (optional)
- ⑥ Vent
- ⑦ Removal cord (optional)
- ⑧

Each hearing instrument can be identified by the serial number located on the outside of the instrument.

* ask your Hearing Care Professional



CIC





(battery type: 10)



Intended use

Hearing instruments are intended to improve the hearing of hearing impaired persons. Diagnosis and prescription of a hearing instrument must be performed by hearing health specialists, e.g. ENT doctors, audiologists or acousticians.

Read and follow the instructions of this user guide to avoid injuries or damage to the hearing instruments.

Symbols:	
Be aware to information marked with the warning symbol WARNING , CAUTION or NOTICE !	
	WARNING points out a situation that could lead to serious injuries, CAUTION indicates a situation that could lead to minor and moderate injuries.
NOTICE	NOTICE indicates possible property damage.
	Advice and tips on how to handle your device better.
	Instruction. Indicates that something has to be done.
	Indicates an option that not every hearing instrument has. Ask your Hearing Care Professional to check off the relevant information.

General notes on safety

NOTICE

Your hearing instruments are sensitive to extreme heat, high humidity, strong magnetic fields ($> 0.1\text{T}$), X-rays and mechanical stress.

- ▶ Do not expose your hearing instruments to extreme temperature or high humidity.
- ▶ Do not leave them in direct sunlight.
- ▶ Do not wear them in the shower or when you apply make-up, perfume, aftershave, hairspray or suntan lotion.
- ▶ Do not wear your hearing instruments when you are exposed to short-waves, a strong magnetic field, a high frequency field or X-rays.
- ▶ Do not place your hearing instruments in a microwave oven.



WARNING

Choking hazard posed by small parts.

- ▶ Keep hearing instruments, batteries and accessories out of children's reach.
- ▶ If swallowed consult a physician or hospital immediately.
- ▶ If infants, small children or mentally disabled persons need to wear hearing instruments ensure adequate supervision.

NOTICE

Leaking batteries damage the hearing instruments.

- ▶ Turn the hearing instruments off when not in use to preserve the battery.
- ▶ Remove batteries when the instruments are not in use for a prolonged period of time.

NOTICE

Batteries contain harmful substances that pollute the environment.

- ▶ Do not throw used batteries into household trash.
- ▶ Dispose batteries according to national regulations or return them to your Hearing Care Professional.



In some countries restrictions for the usage of wireless equipment exist.

- ▶ Refer to local authorities for further information.




WARNING

Hazard of explosion!

- ▶ Do not use your hearing instruments in areas where there is a danger of explosions (e.g. mining).

Handling batteries



original size



- ▶ Always use the right size batteries.
- ▶ Remove the tab on your batteries only when you are ready to use them.
- ▶ The positive (+) side of the batteries is always the smooth side.
- ▶ Always keep spare batteries with you.



- In some types of instruments the sounds may become fainter or you may hear additional sounds. This means that the battery power is low.
- ▶ If this happens, change the battery.

Low battery indicator

Your Hearing Care Professional can activate an acoustic signal to inform you when the battery is nearly empty. If you hear an acoustic signal or experience low amplification power, replace the battery.

Inserting the battery

- ▶ Open the battery compartment.



- ▶ Insert the battery.
Make sure the "+" symbols on the battery and compartment are aligned.



- ▶ Close the battery compartment.

Inserting and removing the hearing instrument



Never pull on the battery compartment door or volume control wheel to remove the hearing instrument, as this could cause damage.


Inserting

- ▶ Insert the hearing instrument into the ear canal.
- ▶ Slightly twist and push the hearing instrument into its proper position.



Removing

- ▶ Push lightly on the back of your ear to loosen the hearing instrument.
- ▶ Pull the hearing instrument out.



Turning your hearing instrument on and off

- ▶ Ask your Hearing Care Professional to check off the relevant information.




Via battery compartment

- ▶ Fully open the battery compartment to turn your hearing instrument off.
- ▶ Fully close the compartment to turn your hearing instrument on.



Via program button

- ▶ Press and hold the button for 3 seconds to turn the instrument on or off.



Via volume control with on/off-function

- ▶ Rotate the volume control down to the stop to turn your hearing instrument off.
- ▶ Rotate the volume control up to turn your hearing instrument on.



Via remote control

- ▶ Refer to the remote control user guide.

Power-on delay

Your hearing instrument is equipped with a "power-on delay" feature which can be activated by your Hearing Care Professional. This allows you to turn the hearing instrument on with an automatic delay of 6, 12 or 18 seconds, so that you may insert the hearing instrument into your ear without experiencing unpleasant feedback (whistling).

Changing the hearing program

- ▶ Ask your Hearing Care Professional to check off the relevant information.



Via program button

- ▶ Press the button to change the program.



Via remote control

- ▶ Refer to the remote control user guide.

Changing the volume

The volume of your hearing instruments is controlled automatically and there is no need to manually adjust it.

Should you prefer manual volume control, your Hearing Care Professional can program the hearing instrument for manual adjustment.




Automatic

The volume is adjusted automatically when in use.



Via volume control

- ▶ Rotate the wheel forward for louder.
- ▶ Rotate the wheel back for quieter.



Via remote control

- ▶ Refer to the remote control user guide.

Functionalities

Your Hearing Care Professional can activate the following functionalities.



e2e wireless

- Allows the use of a remote control.
- Automatically adjusts both hearing instruments simultaneously to change the instrument settings, e.g. volume or program.



The performance of the e2e wireless function may be affected by electromagnetic interference – e.g. from a computer monitor or fluorescent lamp system.

- ▶ Move away from the source of interference if you experience difficulty.



Automatic telephone program

- ▶ Hold a telephone near your hearing instrument (*).

The hearing instrument automatically switches to the telephone program.

- ▶ Move the telephone away from your hearing instrument.

The hearing instrument switches back to the previously used program.



- ▶ You can also select the telephone program manually by using the program button.
- ▶ When operating the hearing instrument in telecoil mode, the wireless function may cause a pulsing noise. This is due to the pick-up of the transmission signal. When your hearing instrument is in telecoil mode please use the remote control only with a distance greater than 10 cm (4 inches).

(*) AutoPhone magnet accessory

Not all telephones produce a magnetic field strong enough to activate the AutoPhone function. In this case, use the Siemens AutoPhone magnet, which is available as accessory, and attach it to the telephone receiver to activate the AutoPhone function.



- ▶ See the AutoPhone magnet user guide for instructions on how to position the magnet.
- ▶ Use only the Siemens magnet. Ask your Hearing Care Professional for this accessory.

NOTICE

Magnets can disturb electrical devices and delete stored data.

- ▶ Keep magnets away from computers, color monitors, television sets, floppy discs, video- and audio discs and other electronic equipment/ devices.



WARNING

Magnets may interfere with the operation of life supporting systems (e. g. pacemakers, magnetic valves).

- ▶ Use a magnet only if not in doubt about the safe distance. The safe distance between pacemaker and magnet should be at least 10 cm (4 inches).

Maintenance and care

Your hearing instruments are highly sensitive systems. Proper care and handling will ensure the best possible performance.

Clean instruments regularly to prevent damage to the hearing instruments and to avoid health problems.

NOTICE

- ▶ Do not put the hearing instruments in water.



- ▶ Clean your hearing instruments with a soft, dry tissue, as necessary.
- ▶ Do not use alcohol or other solvents as it will damage the instruments.
- ▶ Be careful not to damage the wax protection when cleaning the hearing instruments.
- ▶ Place your hearing instruments in a "dry aid kit" with a special moisture-attracting capsule overnight.
- ▶ If the windscreen of your hearing instruments is soiled, ask your Hearing Care Professional to clean or exchange the windscreen.
- ▶ Ask your Hearing Care Professional for special care sets and further information on how to keep your hearing instruments in good condition.

Wax protection

Your hearing instruments are equipped with a sieve, recess or guard on the sound outlet to prevent earwax from entering the hearing instruments.

Your Hearing Care Professional will provide you with the wax protection appropriate for your instruments and can advise you how to clean or replace it.

- ▶ If the instruments sound distorted or weak exchange the wax protection.



Accessories

Tek wireless enhancement

(not all hearing instruments)

Your hearing instruments together with the Tek Connect remote control will allow you to connect your hearing instruments to the growing world of technology. Via a special Bluetooth connection it is possible to receive audio signals from your mobile, your television or your mp3-player.



For detailed information and instructions refer to the Tek user guide or ask your Hearing Care Professional.

Remote control

(not all hearing instruments)

The wireless function of your hearing instruments allows you to use a remote control to change the program and volume of your hearing instruments.

For further information about the remote control, contact your Hearing Care Professional or refer to the remote control user guide.



- ▶ In the event your hearing instruments can be controlled by another person's remote control, return hearing instruments and remote control to your Hearing Care Professional.

Troubleshooting

Use these tips to resolve simple issues.

Problem	Possible solution
Sound is weak.	Replace empty battery. Increase volume. Clean instrument. Close battery door completely or press and hold the program button for approximately 2 seconds.
Hearing instrument emits whistling sounds.	Reinsert the instrument until it fits securely. Clean instrument.
Sound is distorted.	Replace empty battery. Lower volume. Clean instrument.
Hearing instrument emits signal tones.	Gently close battery compartment completely. Replace empty battery.
Hearing instrument does not work.	Check if the hearing instrument is turned on. Gently close battery compartment completely. Replace empty battery. Make sure battery is inserted correctly. Power-on delay is active – wait 18 seconds and check again.

Consult your Hearing Care Professional if you encounter further problems.

Technical information

Siemens Motion 700 WL, Motion 500 WL, Motion 300 WL
Siemens Motion 701 WL, Motion 501 WL, Motion 301 WL
FCC ID: SGI-WL002ITE
IC: 267AB-WL002
Operating frequencies: $F_c=3.28$ MHz



N14203

This Class B digital apparatus complies with Canadian ICES-003.

Changes or modifications not expressly approved by SAT GmbH for compliance could void the users authority to operate the equipment.

FCC compliance

This device complies with Part 15 of the FCC Rules and with RSS-310 of Industry Canada.

Operation is subject to the following conditions:

- this device may not cause harmful interference, and
- this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and

can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Disposal information

Within the European Union, the marked equipment is covered by "Directive 2002/96/EC of the European Parliament and of the Council of 27 January 2003 on waste electrical and electronic equipment."



Amended by "Directive 2003/108/EC" (WEEE).



- ▶ Recycle hearing instruments, batteries and accessories according to national regulations.

With the CE marking Siemens confirms compliance with the European Directive 93/42/EEC concerning medical devices.



For the products with e2e wireless 2.0 Siemens additionally confirms compliance with the European Directive 99/5/EC (R&TTE) concerning radio and telecommunications terminal equipment.

LEFT	Hearing instrument	RIGHT
	Model	
	Serial number	
	Battery	
<input type="checkbox"/> yes <input type="checkbox"/> no	Power-on delay	<input type="checkbox"/> yes <input type="checkbox"/> no
<input type="checkbox"/> yes <input type="checkbox"/> no	Volume control	<input type="checkbox"/> yes <input type="checkbox"/> no
<input type="checkbox"/> yes <input type="checkbox"/> no	Push button	<input type="checkbox"/> yes <input type="checkbox"/> no
Signal tones		
<input type="checkbox"/> on <input type="checkbox"/> off	Low battery warning	<input type="checkbox"/> on <input type="checkbox"/> off
<input type="checkbox"/> on <input type="checkbox"/> off	Volume control limit	<input type="checkbox"/> on <input type="checkbox"/> off
<input type="checkbox"/> on <input type="checkbox"/> off	volume control change	<input type="checkbox"/> on <input type="checkbox"/> off
Hearing program	Description of the hearing program	
1		
2		
3		
4		
5		
Remote control		
Model		
Serial number		

© Siemens AG, 04.2010 · SB
Document No. A91SAT-00673-99T5-7600
Order/Item No. 104 202 05
Printed in Germany

Siemens Audiologische Technik GmbH
Gebbertstrasse 125
91058 Erlangen
Germany
Phone +49 9131 308 0

www.siemens.com/hearing